

SiliBand®

Gastroplastie Verticale Calibrée

Anneau en silicone

- > Calibrage de la poche de gastroplastie
- > Minilaparotomie.

Ring made of silicone

- > *Calibration of the gastroplasty pouch*
- > *Minilaparotomy*




**SURGICAL
IOC®**

L'alliance pour l'innovation

SiliBand®

Gastroplastie Verticale Calibrée




Positionnement de l'anneau sur l'estomac
Positioning of the ring on the stomach



Fixation de l'anneau
Fixing of the ring to the stomach



Cicatrice < 7 cm
Scar < 7 cm



Contrôle radio post-opératoire
Post-operative radiography

Anneau Siliband® :

- Constitué de Silicone non-adhérent et radio-opaque.
- Facile à positionner sur l'estomac.
- Forme et géométrie garantissant un calibrage stable et reproductible de la poche gastrique.

Composition du kit :

- 1 anneau SiliBand® pour Gastroplastie en silicone, opaque aux rayons X
- 1 fil de suture 2/0 non résorbable tressé en polyester revêtu anti-adhérent avec aiguille droite de 60 mm, pointe ronde
- 1 pince de préhension
- 1 seringue de 10cc avec rince cathéter
- 1 règle centimétrique

Documents annexes :

- 1 notice d'utilisation
- 1 vignette pour identification et traçabilité
- 1 livret patient

Réf.	L (mm)	Ø* (mm)
SB 500	50	15,7
SB 550	55	17,5

*diamètre interne de l'anneau fermé
*internal diameter of the closed ring

Siliband® the Ring :


- Made of non-adherent and radio-opaque Silicone.
- Easy to set on the stomach.
- Shape and geometry ensuring a stable and reproducible calibration of the gastric pouch.

Composition of kit :

- 1 SiliBand® Gastroplasty ring made of X-ray opaque silicone
- 1 non-absorbable interlaced 2/0 suture thread made of anti-adherent coated polyester, with straight, rounded tip, 60 mm needle
- 1 grasping tweezer
- 1 10cc syringe
- 1 centimetric ruler

Additional documents :

- 1 instructions for use
- 1 sticker for identification and traceability
- 1 patient leaflet



Bibliographie :

> Pierre Urbain, MD. The Proring®-Band : A New Device in the Field of Restrictive surgery : Preliminary Experience. *Obes Surg* 2002 ; 12 : 588-591

> Commandes | Orders :

Tél. +33(0)4 77 91 42 00 | Fax. +33(0)4 77 91 42 01

SURGICAL - IOC® - 88 rue Marengo
42000 SAINT - ETIENNE | FRANCE
Internet : www.surgical-ioc.com



**SURGICAL
IOC®**
L'alliance pour l'innovation